

Bedside cabinets



Cherusker

- pivoting overbed table
- height adjustment of the overbed table 82–112 cm
- width 52 cm, overall width (table horizontal) 114 cm
- double castors Ø 50 mm, 2 locking castors



Arminius

- height-adjustable and pivoting overbed table
- width 54,5 cm, overall width (table horizontal) 125 cm
- double castors Ø 50 mm, 2 locking castors



Hermann

- retractable overbed table
- width 55 cm, overall width (table horizontal) 119,5 cm
- double castors Ø 50 mm, 2 locking castors



Hermann

- without overbed table
- width 44,5 cm
- double castors Ø 50 mm, 2 locking castors



Varus

- width 54 cm, height 57 cm
- double castors Ø 50 mm
- available in Natural Beech, Royal Maple, Cherry Havana and Bella Noce Choco
- ideally suited for low-height beds

Servers, Ravello, Relax



Duo

- pivoting overbed table
- height adjustment of the overbed table 80–113 cm
- fixed side tray
- double castors Ø 50 mm, 2 locking castors



Solo

- pivoting overbed table
- height adjustment of the overbed table 80–113 cm
- double castors Ø 50 mm, 2 locking castors



Ravello

- model variant Curo: fixed position footrest tilts to follow backrest adjustments
- model variant Sano: footrest can be fully retracted
- backrest inclines to a 43° reclining position



Relax

- homelike wooden bed surround
- for mattress bases 90 x 200, 100 x 200, 120 x 200, 120 x 220, 140 x 200, 140 x 220 cm
- combinable with motorised frame and the Lippe bed inserts
- head and footboards with frames made of solid beech wood, 4 décors available for fillings and side panels



▲ Bella Noce Choco



▲ Natural Beech



▲ Cherry Havana



▲ Royal Maple

Design



Softcovers

- textile coatings in 7 attractive colours for headboards, footboards and safety sides
- easy to apply, fully launderable
- for care beds of the Dali range (Dali low-econ, Dali low-entry-econ, Dali low-entry: only safety sides)

▼ Your retailer

For all questions regarding care beds, please consult your competent specialist dealer



BURMEIER

Burmeier GmbH & Co. KG
Pivitsheider Straße 270
32791 Lage
Phone +49 (0) 5232 9841 - 0
Fax +49 (0) 5232 9841 - 41
info@burmeier.com
www.burmeier.com

A Stieglmeyer-Group company

STIEGLMEYERGROUP

burmeier.com

SAFE LIVING AT HOME



Beds and furniture for care at home



International edition

Subject to technical changes without notice. The shown product colours may differ slightly from the actual product. Last updated 03/2019, Version 6, Item no. 600993.

STIEGLMEYERGROUP



BURMEIER



Burmeier – safe living at home

‘Safe living at home’ – this is the philosophy that characterises Burmeier products. The company, with its roots in Lage in the Westphalian region of Lippe, is the market leader in the home care sector in Germany. The technically innovative and durable products enjoy a worldwide reputation.

Burmeier has been another successful subsidiary of the Stiegmeyer Group in its current form since 1995. Product development always focuses on the concept of making life at home as pleasant as possible for people in need of care and making work easier for care-giving relatives.

Attractive design, modern safety side systems and prompt availability of the products are key arguments for retailers and private customers. The Burmeier 24-volt drive system does not require any further leakage current measurements during the expected service life of the bed and beyond. All this saves costs and helps to ensure that the products offer excellent value for money.

In Burmeier’s safe, comfortable beds you or your customers will feel very comfortable.

Care beds



Dali econ

- mattress base width 90 cm
- safe working load 185 kg
- height adjustment 40–80 cm
- safety sides with Easy Click
- optional wooden encasements
- Econ 230-volt drive system



Dali low-econ

- safe working load 185 kg
- height adjustment 32–72 cm
- optional wooden encasements
- Econ 230-volt drive system

Dali low-entry-econ

- safe working load 175 kg
- height adjustment 23–63 cm
- optional wooden encasements
- Econ 230-volt drive system



Dali

- mattress base width 90 cm
- safe working load 185 kg
- height adjustment 40–80 cm
- wireless handset
- safety sides with Easy Click
- optional wooden encasements
- 24-volt drive system
- Dali-Lock app



Dali wash

- mattress base width 90 cm
- safe working load 185 kg
- height adjustment 40–80 cm
- automatic reprocessing
- wireless handset
- safety sides with Easy Click
- 24-volt system
- Dali-Lock app



Dali low-entry

- mattress base width 90 cm
- safe working load 175 kg
- height adjustment 23–63 cm
- wireless handset
- safety sides with Easy Click
- optional wooden encasements
- 24-volt drive system
- Dali-Lock app

Residence and senior beds



Westfalia IV

- mattress base width 80, 90, 100 cm
- safe working load 185 kg
- height adjustment 35–80 cm
- several wooden surrounds to choose from
- 24-volt drive system



Westfalia-Klassik

- mattress base width 90 cm
- safe working load 200 kg
- height adjustment 35–80 cm
- central braking system
- several wooden surrounds to choose from
- 24-volt drive system



Westfalia-Klassik Elegant variant

- light Magnolia décor
- Lito head and footboard with homelike curved shape
- framework colour Argentum
- concealed castors



Easy
Switch

Regia

- mattress base width 90, 100 cm
- safe working load 225 kg
- height adjustment 26–80 cm
- telescopic safety sides can be raised in 2 stages
- several wooden surrounds to choose from
- 24-volt drive system
- modular system Easy Switch



Inovia

- mattress base width 90, 100 cm
- safe working load 225 kg
- height adjustment 22–77 cm
- full-length safety sides
- telescopic safety sides optional
- several wooden surrounds to choose from
- 24-volt drive system

Mobilisation and special beds



Lippe IV + Relax

- mattress base width 80, 90, 100 cm (combinable with bed surround Relax: 90, 100 cm)
- safe working load 200 kg
- height adjustment 29–74 cm
- 24-volt drive system
- machine-washable variant optional



Lippe 120/140 + Relax

- mattress base width 120 and 140 cm (both combinable with bed surround Relax)
- safe working load 250 kg
- height adjustment 29–74 cm
- 24-volt drive system



Allura II

- mattress base width 100, 120 cm
- safe working load 250 kg
- height adjustment 30–80 cm
- reverse Trendelenburg position electrically adjustable
- 24-volt drive system



Gigant

- mattress base width 120 cm
- safe working load 350 kg
- height adjustment 40–80 cm
- reverse Trendelenburg position electrically adjustable
- split safety sides optional



Vertica homecare

- mattress base width 90 cm
- safe working load 225 kg
- height adjustment 40–80 cm
- moves the resident in an upright sitting position at the touch of a button