



The heavy-duty care bed


Gigant

ADVANTAGES

The Gigant heavy-duty bed can support occupants weighing up to 280 kg, making it suitable for extremely obese persons. It also offers plenty of space and support for independent living and reduces the burden on carers.

- Safe working load: 350 kg
- Maximum weight of resident: 280 kg
- Height adjustment range: 40 to 80 cm
- Reverse-Trendelenburg position
- Optional combinable safety side
- Headboards and footboard with elegant wooden surround

External dimensions/height adjustment range:



Special features:



Further information:



DESIGN DETAILS



Holder for mobilisation aid (left) and locking box (right).



CPR release of the backrest.



Extremely sturdy safety side construction.



Plastic panels on the chassis for a more homelike feel.

ACCESSORIES



Clamp-on safety side bar.



Combinable safety side.



Reading lamps.



Bed extension including longer safety sides and side panels.



Foam leather cover for telescopic safety sides (TSG).



Extension cushion to place over the full-length safety side.



The storage tray for glasses, mobiles or the wireless handset can be attached at various positions.

EQUIPMENT AND OPTIONS

Equipment/Option	Gigant	Material no.
Design features		
Gigant 120 x 200 cm		238647LFM013
Technical data		
Net weight	159 kg	
Overall dimensions	135 x 222 cm	
Height adjustment range	40 – 80 cm	
Safe working load	350 kg	
Maximum resident weight	280 kg	
Standard features		
Frame colour	Argentum	
Decor	Natural Beech	
Mattress base size	120 x 200 cm	
Safety side	Full-length safety side (DSG)	
Type of handset (selectively lockable / wireless)	Selective	
Mechanical CPR release	✓	
Castor diameter	100 mm	
Castor locking	Single	
Ratchet mechanism, electrical adjustment of lower leg rest	✓	
Accessories		
Clamp-on safety side bar	O	164977
Combinable safety side on both sides	O	SEG022
Combinable safety side on left side	O	SEG023
Combinable safety side on right side	O	SEG024
Bed extension for metal slatted mattress base (with DSG)	O	240271 + 242308
Mattress piece for bed extension 120 cm wide, 19 cm thick	O	712167
Foam leather cover for DSG, one side	O	225921
Foam leather cover for KSG, one side	O	2x 183008 + 1x 712961
Reading lamp for easy insertion	O	264877
Reading lamp for fixing to the patient lifting pole	O	1x 191879 1x 174231
Extension cushion (for full-length safety side only)	O	258395
Storage tray	White Yellow Red	266970 266971 266619