



Heavy duty beds

Room to stretch out

Gigant



Allura II



Lippe 120/140





Strong

“In my Allura II, with its 120 cm-wide mattress base, I wake up on a morning more rested and refreshed than ever before. Now I’m finally able to start the day full of vitality again.”

Agnes Brown
Customer



Allura II



The cosy heavy-duty bed

For decades, people have been growing taller and heavier, on average. This trend also places greater demands on care provision. Care beds with a wider mattress base and a higher load capacity ensure greater comfort and protection. They are also particularly future-proof.

- Safe working load of 250 kg
- Allura II 100 and Allura II 120 models with mattress-base widths of 100 and 120 cm
- Optional bed extension for a mattress-base length of 220 cm
- Optionally available standing aid promotes the resident’s mobility
- Reverse-Trendelenburg position for a comfortable sitting position
- Full-length Easy Click safety sides with tool-free bar attachment for effortless installation and thorough cleaning

Gigant



The bed for the most demanding requirements

The Gigant heavy-duty bed offers maximum comfort and support even for people living with marked obesity. Its height adjustment range of 40 to 80 cm also makes bedside care much easier. The optional combinable safety sides provide convenient support for the resident when getting into and out of bed.

- Safe working load of 350 kg
- Suitable for residents weighing up to 280 kg
- 120 cm-wide mattress base
- Optional bed extension for a mattress base length of 220 cm
- Optional combinable safety sides with mobilisation aid for protection and greater independence
- Headboards and footboards with elegant wooden surround
- Side panels for a cosy ambience



Twosome

“By combining the Relax with the Lippe, we have created our ideal double bed – and it fits right in with our furnishings.”

*Jacob and
Olivia Jones*
Married couple



Lippe 120/140



Height-adjustable bed inserts

The Lippe 120 and 140 bed inserts can be easily installed in an existing bed surround to provide the comfort and safety of a first-class care bed.

- Safe working load of 250 kg
- 120 or 140 cm wide mattress base
- Optional bed extension for a mattress base length of 220 cm
- Optional standing aid promotes the resident's mobility
- Combines perfectly with our stylish Relax bed surround

Nicht alle Produkte und Konfigurationen sind in jedem Land verfügbar. Not all products and configurations are available in every country. Technische Änderungen vorbehalten. Leichte Farbabweichungen der Originale sind drucktechnisch bedingt. Stand 12/2024. Version 04-125. Artikel-Nr. 601745. Subject to technical changes without notice. The shown product colours may differ slightly from the actual product. Last updated 12/2024. Version 04-125. Item no. 601745.



stieglmeyer-group.com

Burmeier GmbH & Co. KG
Industriestraße 53
32120 Hiddenhausen

Phone +49 (0) 5223 9769 - 0
Fax +49 (0) 5223 9769 - 090
info@burmeier.com