

## SAFE LIVING AT HOME



The care bed for a more homelike ambiance


# Westfalia IV

## ADVANTAGES


The Westfalia IV combines excellent comfort, advanced safety and a particularly stylish design. The height can be adjusted with the help of a stable yet sleek and elegant scissor mechanism. The headboards and footboards are available in various attractive forms.

- Safe working load: 185 kg
- Maximum weight of resident: 155 kg
- Height adjustment range: 35 to 80 cm
- 24-volt drive system
- No leakage current measurements during expected service life
- Selectively lockable handset
- Attractive head and footboard designs
- Different mattress base sizes to choose from

External dimensions/height adjustment range:



Special features:



### Further information:



Product page



Assembly video

## The advantages at a glance

## DESIGN DETAILS



Handset with selective locking function.



Individually braked castors.



Height adjustment range 35 to 80 cm.



Patient lifting pole.



Choice of headboards and footboards: Standard (standard feature), Comodo with round posts (optional), Primero with handle bar (optional).

## ACCESSORIES



Tray.



Mobilisation aid.



Adaptable safety sides.



Wall deflection roller including mount.



Reading lamps.



Bed extension including safety side bars (also for retrofitting).



Foam leather cover.



The storage tray for glasses, mobiles or the wireless handset can be attached at various positions.



Flexible transmitter holder for the handset.



Extension cable for the handset (1.40 m).



With the floor-mounted power supply unit, the 24-volt system can also be used with various international mains cables.

## EQUIPMENT AND OPTIONS

Equipment/Option	Westfalia IV	Material no.
<b>Design features</b>		
Westfalia IV 90 x 200 cm		229911ZLFB000
Westfalia IV 100 x 200 cm <sup>1)</sup>		229911ZLFB002
Westfalia IV 90 x 180 cm		229911ZLFL001
Westfalia IV 90 x 190 cm		229911ZLFL010
<b>Technical data</b>		
Net weight	136 kg	
Overall dimensions	110 x 232 110 x 212 100 x 232 100 x 192 100 x 202 100 x 212 90 x 212	
Height adjustment range	35 – 80 cm	
Safe working load	185 kg	
Maximum resident weight	155 kg	
24-volt system	✓	
<b>Standard features</b>		
Frame colour	Topaz	
Decor	Natural Beech	
Mattress base size	90 x 200 cm	
Sprung wooden slat mattress base (mattress base model)	✓	
Safety side	Full-length safety side (DSG)	
Type of handset (selectively lockable/wireless)	Selective	
Mechanical CPR release	✓	
Castor diameter	125 mm	
Castor locking	Single	
Ratchet mechanism, electrical adjustment of lower leg rest	✓	
Patient lifting pole	✓	
<b>Special equipment and accessories</b>		
Adaptable safety sides	Right Left	250138 250139
Mobilisation aid	O	175402
Comodo wooden surround with round posts	O	ZMOK015
Primero wooden surround	O	ZMOK017
Bed extension to mattress base length of 220 cm	O	ZLFL007
Mattress piece for bed extension 90 cm wide, 12 cm thick	O	166715
Mattress piece for bed extension 100 cm wide, 12 cm thick	O	219857
4 wall deflection rollers, horizontal action	O	ZWAN004
2 wall deflection rollers, horizontal action, foot end	O	ZWAN003
2 wall deflection rollers, horizontal action, head end	O	ZWAN002
Foam leather cover for full-length safety sides (DSG)	O	225921
Reading lamp for easy insertion	O	264877
Reading lamp for fixing to the patient lifting pole	O	1 x 191879 1 x 174231
Tray for full-length safety sides (DSG), 90 cm wide	O	174296
Tray for full-length safety sides (DSG), 100 cm wide	O	273465
Extension cushion (for full-length safety side only)	O	258395
Transmitter holder for handset	O	199137
Handset extension cable	O	178986
Document holder	O	204241
Storage tray	White Yellow Red	266970 266971 266619
Floor-mounted 24-volt power supply unit	O	254276

✓ = Standard feature O = Available as an option <sup>1)</sup> With Comodo head/footboard only