



The care bed
with special features


Westfalia-Klassik

ADVANTAGES


The Westfalia-Klassik is designed for discerning users who are not prepared to compromise on comfort, individuality and sophisticated technical features in their daily care routine.

- Safe working load: 200 kg
- Maximum weight of resident: 170 kg
- Height adjustment range: 35 to 80 cm
- 24-volt drive system
- No leakage current measurements during expected service life
- Sliding backrest in accordance with the recommendations of the German Nurses Association (DBfK)
- Central castor locking

External dimensions/height adjustment range:



Special features:



Further information:



Product page

The advantages at a glance

DESIGN DETAILS



Central braking possible whatever the mattress base position.



Electronic (SwitchMode) power supply unit integrated immediately downstream of the mains plug.



Mattress base height adjustable from 35 to 80 cm, adjustment range = 45 cm.

ACCESSORIES



Tray.




Wall deflection roller.



Wooden surrounds: Standard (standard feature), Comodo with round posts (optional), Primero with handle bar (optional).



Reading lamps.



Bed extension with metal slats, safety sides and side panels.



Foam leather cover.



The storage tray for glasses, mobiles or the wireless handset can be attached at various positions.



Flexible transmitter holder for the handset.



Extension cable for the handset (1.40 m).



Extension cushion to place over the full-length safety side.



With the floor-mounted power supply unit, the 24-volt system can also be used with various international mains cables.

EQUIPMENT AND OPTIONS

Equipment/Option	Westfalia-Klassik	Material no.
Design features		
Westfalia-Klassik 90x200 cm with standard wooden surround		171358HOL006
Technical data		
Net weight	115 kg	
External dimensions with standard wooden surround	102 x 208 cm	
Height adjustment range	35 – 80 cm	
Safe working load	200 kg	
Maximum resident weight	170 kg	
24-volt system	✓	
Standard features		
Frame colour	Topaz	
Decor	Natural Beech	
Mattress base size	90 x 200 cm	
Safety side	Full-length safety side (DSG)	
Type of handset (selectively lockable / wireless)	Selective	
Mechanical CPR release	✓	
Castor diameter	100 mm	
Castor locking	Central	
Ratchet mechanism, electrical adjustment of lower leg rest	✓	
Accessories		
Comodo wooden surround with round posts	O	HOL007
Primero wooden surround	O	HOL012
Bed extension	O	701215 + 203806
Wall deflection roller Ø 100 mm	O	WAN011
Foam leather cover for full-length safety sides (DSG)	O	225921
Reading lamp for easy insertion	O	264877
Reading lamp for fixing to the patient lifting pole	O	1 x 191879 1 x 174231
Tray for full-length safety sides (DSG)	O	174296
Extension cushion (for full-length safety side only)	O	258395
Transmitter holder for handset	O	199137
Handset extension cable	O	178986
Document holder	O	204241
Storage tray	White	266970
	Yellow	266971
	Red	266619
Floor-mounted 24-volt power supply unit	O	204276