

## Bedside cabinets



### Cherusker

- pivoting overbed table
- height adjustment of the overbed table 82–112 cm
- width 52 cm, overall width (table horizontal) 114 cm
- double castors Ø 50 mm, 2 locking castors



### Arminius

- height-adjustable and pivoting overbed table
- width 54,5 cm, overall width (table horizontal) 125 cm
- double castors Ø 50 mm, 2 locking castors



### Hermann

- retractable overbed table
- width 52,5 cm, overall width (table horizontal) 114 cm
- double castors Ø 50 mm, 2 locking castors



### Hermann

- without overbed table
- width 44,5 cm
- double castors Ø 50 mm, 2 locking castors



### Varus

- width 54 cm, height 57 cm
- double castors Ø 50 mm
- available in Natural Beech, Royal Maple, Havana Cherry and Bella Noce Choco
- ideally suited for low-height beds

## Servers, Ravello, Relax



### Duo

- pivoting overbed table
- height adjustment of the overbed table 80–113 cm
- fixed side tray
- double castors Ø 50 mm, 2 locking castors



### Solo

- pivoting overbed table
- height adjustment of the overbed table 80–113 cm
- double castors Ø 50 mm, 2 locking castors



### Ravello

- model variant Curo: fixed position footrest tilts to follow backrest adjustments
- model variant Sano: footrest can be fully retracted
- backrest inclines to a 43° reclining position



### Relax

- homelike wooden bed surround
- for mattress bases 90 x 200, 100 x 200, 120 x 200, 120 x 220, 140 x 200, 140 x 220 cm
- combinable with motorised frame and the Lippe bed inserts
- head and footboards with frames made of solid beech wood, 4 décors available for fillings and side panels



▲ Bella Noce Choco



▲ Natural Beech



▲ Havana Cherry



▲ Royal Maple

## Design



### Softcovers

- textile coatings in 6 attractive colours for headboards, footboards and safety sides
- easy to apply, fully launderable
- for care beds of the Dali range (Dali low-econ, Dali low-entry-econ, Dali low-entry: only safety sides) and for our Lenus (safety sides)



### Comfort-fit care bed mattress



### Mattresses

For restful sleep, the choice of the mattress is just as important as the bed. Our mattresses optimally fulfil the various needs of homecare. Learn more about our mattresses by scanning the QR code or directly at [www.burmeier.com/en/mattresses](http://www.burmeier.com/en/mattresses).



▼ Your retailer

For all questions regarding care beds, please consult your competent specialist dealer



Burmeier GmbH & Co. KG  
Industriestraße 53  
32120 Hiddenhausen  
Phone +49 (0) 5223 9769 - 0 Fax  
+49 (0) 5223 9769 - 090  
info@burmeier.com  
www.burmeier.com

A Stieglmeier-Group company

STIEGLMEYERGROUP

burmeier.com

## SAFE LIVING AT HOME

Subject to technical changes without notice. The shown product colours may differ slightly from the actual product. Last updated 05/2021. Version 07. Item no. 601303.



*Beds and furniture for care at home*



International edition

STIEGLMEYERGROUP





## Burmeier – safe living at home

‘Safe living at home’ – this is the philosophy that characterises Burmeier products. The company, with its roots in Hiddenhausen in Westphalia, is the market leader in the homecare sector in Germany. The technically innovative and durable products enjoy a worldwide reputation.

Burmeier has been another successful subsidiary of the Stieglmeyer Group in its current form since 1995. Product development always focuses on the concept of making life at home as pleasant as possible for people in need of care and making work easier for care-giving relatives.

Attractive design, modern safety side systems and prompt availability of the products are key arguments for retailers and private customers. The Burmeier 24-volt drive system does not require any further leakage current measurements during the expected service life of the bed and beyond. All this saves costs and helps to ensure that the products offer excellent value for money.

In Burmeier’s safe, comfortable beds you or your customers will feel very comfortable.

### Care beds



#### Dali econ

- mattress base width 90 cm
- safe working load 185 kg
- height adjustment 40–80 cm
- safety sides with Easy Click
- optional wooden encasements
- Econ 230-volt drive system



#### Dali low-econ

- safe working load 185 kg
- height adjustment 32–72 cm
- optional wooden encasements
- Econ 230-volt drive system



#### Dali

- mattress base width 90 cm
- safe working load 185 kg
- height adjustment 40–80 cm
- wireless handset
- safety sides with Easy Click
- optional wooden encasements
- 24-volt drive system



#### Dali wash

- mattress base width 90 cm
- safe working load 185 kg
- height adjustment 40–80 cm
- automatic reprocessing
- wireless handset
- safety sides with Easy Click
- 24-volt system



#### Dali low-entry

- mattress base width 90 cm
- safe working load 175 kg
- height adjustment 23-63 cm
- wireless handset
- safety sides with Easy Click
- optional wooden encasements
- 24-volt drive system

### Residence and senior beds



#### Westfalia IV

- mattress base width 80, 90, 100 cm
- safe working load 185 kg
- height adjustment 35–80 cm
- several wooden surrounds to choose from
- 24-volt drive system



#### Westfalia-Klassik

- mattress base width 90 cm
- safe working load 200 kg
- height adjustment 35–80 cm
- central braking system
- several wooden surrounds to choose from
- 24-volt drive system



#### Inovia

- mattress base width 90, 100 cm
- safe working load 225 kg
- height adjustment 22–77 cm
- full-length safety sides
- telescopic safety sides optional
- several wooden surrounds to choose from
- 24-volt drive system



#### Regia

- mattress base width 90, 100 cm
- safe working load 225 kg
- height adjustment 26–80 cm
- telescopic safety sides can be raised in 2 stages
- several wooden surrounds to choose from
- 24-volt drive system
- modular system Easy Switch



#### Regia partner

- mattress base width 2 x 90 cm
- safe working load 2 x 225 kg
- height adjustment 26 - 80 cm
- separable double bed
- modular system Easy Switch
- comfortable sitting position
- wooden surrounds and homelike upholstery to choose from

### Special beds



#### Lippe IV + Relax

- mattress base width 80, 90, 100 cm (combinable with bed surround Relax: 90, 100 cm)
- safe working load 200 kg
- height adjustment 29–74 cm
- 24-volt drive system
- machine-washable variant optional



#### Lippe 120/140 + Relax

- mattress base width 120 and 140 cm (both combinable with bed surround Relax)
- safe working load 250 kg
- height adjustment 29–74 cm
- 24-volt drive system



#### Allura II

- mattress base width 100, 120 cm
- safe working load 250 kg
- height adjustment 30–80 cm
- reverse Trendelenburg position electrically adjustable
- 24-volt drive system



#### Gigant

- mattress base width 120 cm
- safe working load 350 kg
- height adjustment 40–80 cm
- reverse Trendelenburg position electrically adjustable
- split safety sides optional



#### Lenus

- mattress base width 90 cm
- safe working load 225 kg
- height adjustment 15 - 80 cm
- safety sides with Easy Click
- optional wooden encasements
- 24-volt drive system