



Practical and functional

Bedside cabinets

Further information:



Product page

Bedside cabinets and side tables

Cherusker bedside cabinet

With height-adjustable tilting table, plastic-coated body, beech decor

Technical data:

Width: 52 cm
Height: 87 cm
Overall depth: 50 cm
Overall width (with table horizontal): 114 cm
Table height adjustment from 82 to 112 cm
Double castors Ø 50 mm, 2 castors are lockable



Item number: 193753



Arminius bedside cabinet

With height-adjustable tilting table, plastic-coated body, beech decor

Technical data:

Width: 54.5 cm
Height: 89 cm
Overall depth: 45 cm
Overall width (with table horizontal): 125 cm
Table height adjustment from 78 to 112 cm
Double castors Ø 50 mm, 2 castors are lockable



Item number: 236959





Item number: 223636

Hermann bedside cabinet with table

With fold-down table, plastic-coated body, beech decor

Technical data:

Width: 55 cm
Height: 89 cm
Overall depth: 45 cm
Overall width (with table horizontal): 119.5 cm
Double castors Ø 50 mm, 2 castors are lockable



Item number: 232644

Varus bedside cabinet

Plastic-coated body, no table, available in four different decors: Natural Beech, Royal Maple, Cherry Havana, Bella Noce Choco

Technical data:

Width: 54 cm
Height: 57 cm
Overall depth: 45 cm
Double castors Ø 50 mm, 2 drawers



Item number: 218931

Duo overbed table and side table

With height-adjustable, tilting overbed table and fixed table section, beech decor

Technical data:

Width: 81 cm
Height: 81 cm
Overall depth: 42 cm
Overbed table height adjustable from 80 to 113 cm
Double castors Ø 50 mm, 2 castors are lockable



Item number: 243335

Hermann bedside cabinet without table

Plastic-coated body, beech decor, no table

Technical data:

Width: 44.5 cm
Height: 89 cm
Overall depth: 45 cm
Double castors Ø 50 mm, 2 castors are lockable



Item number: 218932

Solo overbed table and side table

With height-adjustable tilting table, beech decor

Technical data:

Width: 60 cm
Height: 81 cm
Overall depth: 42 cm
Overbed table height adjustable from 80 to 113 cm
Double castors Ø 50 mm, 2 castors are lockable