



The care chair for extra mobility

Ravello

ADVANTAGES

The Ravello transport chair allows persons in need of care to be transferred safely from the bedroom to the bathroom or dining room, for example. Its flexible construction makes it easy for carers to transfer the resident to the chair from an armchair or the side of the bed. The cushioned, ergonomically shaped seat also makes the Ravello comfortable to sit in for longer periods. The washable seat covers are very easy to keep clean.

- Safe working load of 135 kg (dynamic), 200 kg (static)
- The lowerable armrests make it easy to get in to the chair from the side.
- The backrest can be lowered to a 43° reclining position.
- Ravello-sano: retractable footrest
- Ravello-curo: permanently fixed footrest moves with the lower leg rest during adjustments
- Ravello sano-s: 4 steerable castors for easy manoeuvring
- Ravello plus-s: height-adjustable seat

Further information:



The advantages at a glance

DESIGN DETAILS



Ravello sano: the footrest can be completely retracted.



Ravello curo: length-adjustable, fixed position footrest that tilts to follow backrest adjustments.



Armrest height adjustment in 6 stages.



The Ravello sano-s is available for better manoeuvrability.

ACCESSORIES



Side supports.



Headrest with lateral support.



Urine bag holder.



Attachable table.



Infusion stand.

EQUIPMENT AND OPTIONS

Equipment/Option	Ravello	Material no.
Design features		
Ravello curo		264299
Ravello sano		264297
Ravello sano-s		264298
Ravello plus-s		264300
Technical data		
Net weight	Curo, Sano, Sano-s: 38 kg Plus-s: 58 kg	
Dimensions	Curo, Sano, Sano-s: W x H x D 68 x 126 x 103 cm incl. footrests Plus-s: W x H x D 73 x 132-161 x 109 cm incl. footrests	
Backrest	Height 81 cm, reclines to 43°	
Footrest	26 x 24 cm, with 8 cm ground clearance	
Castors	Front castors: 20 cm, fixed Rear castors: 15 cm, swivel type, with paired locking	
S version castors:	15 cm, swivel type, with paired locking	
Seat	PU foam with wooden core	
Backrest	PU foam with steel core	
Covering fabric	Faux leather	
Special equipment and accessories		
Tray, attachable table	0	195900
Neck support (black)	0	116344
1 pair of side supports (azure blue)	0	189619
Leg support, right	0	116342
Leg support, left	0	116343
Urine bag holder	0	116340
Infusion stand	0	1x 260269 1x 114352